

## Sterilisation and Instructions for Use – Product Schedule

**Product: Prophy Cups, RA, Screw and Snap-On**  
**Classification: Class 2a Medical Device – (single use device)**

GMDN Code	GMDN Term	Basic UDI-DI	Medical Device	MD Class
46938	Dental Polishing Cup, single use	5051717PC01UF	Prophy Polishers	Ila

PRA001	Prophy Cup, Pink, 6 Web, Ring, RA
PRA002	Prophy Cup, Purple, 6 Web, Ring, RA
PRA003	Prophy Cup, Blue, 6 Web, Ring, RA
PRA004	Prophy Cup, White, 6 Web, RA
PRA005	Prophy Cup, Yellow, 6 Web, RA
PSRA001	Prophy Cup, Pink, 6 Web, RA
PSRA002	Prophy Cup, Purple, 6 Web, RA
PSRA003	Prophy Cup, Blue, 6 Web, RA
PSRA004	Prophy Cup, White, RA
PSRA005	Prophy Cup, Yellow, RA
PSRA011	Prophy Cup, Pink, 4 Web, Ribbed, RA
PSRA012	Prophy Cup, Purple, 4 Web, Ribbed, RA
PSRA013	Prophy Cup, Blue, 4 Web, Ribbed, RA
PSRA014	Prophy Cup, Blue, 4 Web, Ribbed, RA
PSRA015	Prophy Cup, Yellow, 4 Web, RA
PSRA017	Prophy Cup, Green, 4 Web, RA
TPY002G	Prophy Cup, White, 4 Web, RA
EPR003	Prophy Cup, Mint, 6 Web Ring, RA
EPY002	Prophy Cup, Mint, 4 Web, Ribbed, RA
EPY003A	Prophy Cup, Mint, 6 Web, RA
PSC001	Prophy Cup, Pink, 6 Web Screw
PSC002	Prophy Cup, Purple, 6 Web Screw
PSC003	Prophy Cup, Blue, 6 Web Screw
PSC004	Prophy Cup, White, 6 Web Screw
PSC005	Prophy Cup, Yellow, 6 Web Screw
PSC011	Prophy Cup, Pink, 4 Web, Ribbed Screw
PSC012	Prophy Cup, Purple, 4 Web, Ribbed Screw
PSC013	Prophy Cup, Blue, 4 Web, Ribbed Screw
PSC014	Prophy Cup, Blue, 4 Web, Ribbed Screw
PSC015	Prophy Cup, Yellow, 4 Web, Ribbed Screw
PSC015XH	Prophy Cup, Yellow, 4 Web, Ribbed Screw
PSC018	Prophy Cup, Green, 4 Web, Ribbed, Screw
EPR002A	Prophy Cup, White, 6 Web Screw
PSMRA001	Prophy Cup Mini, Pink, 4 Web, Ribbed, RA
PSMRA002	Prophy Cup Mini, Purple, 4 Web, Ribbed, RA
PSMRA003	Prophy Cup Mini, Blue, 4 Web, Ribbed, RA
PSMSC001	Prophy Cup Mini, Pink, 4 Web, Ribbed, Screw
PSMSC002	Prophy Cup Mini, Purple, 4 Web, Ribbed, Screw
PSMSC003	Prophy Cup Mini, Blue, 4 Web, Ribbed, Screw
EP001	Prophy Cup, White, 4 Web, Ribbed, RA
EPY002G	Prophy Cup, Mint, 4 Web, RA
EPY004	Prophy Cup, Grey, Ribbed, RA
EPY005	Prophy Cup, Grey, Ribbed, Screw
TPR002A	Prophy Cup, White, 6 Web, Screw
TPR002B	Prophy Cup, 4 Web, Screw
TPR002C	Prophy Cup, Blue, 4 Web, Screw
TPR002E	Prophy Cup, Yellow, 4 Web, Ribbed, Screw
TPR002M	Prophy Cup, Grey, 4 Web, Screw
TPR002H	Prophy Cup, Green, 4 Web, Screw
TPR003	Prophy Cup, White, 6 Web, Ring, RA
TPR003B	Prophy Cup, Blue, 4 Web, RA

## Sterilisation and Instructions for Use

TPR003E	Prophy Cup, White, Ribbed, RA
TPR003D	Prophy Cup, Yellow, 4 Web, Ribbed, RA
TPR003G	Prophy Cup, Green, 4 Web, RA
TPR020	Prophy Cup, Yellow, 4 Web, Ribbed, Screw
TPR021	Prophy Cup, Grey, 4 Web, Ribbed, Screw
TPR022	Prophy Cup, Green, 6 Web, Screw
TPR023	Prophy Cup, Yellow, 4 Web, Ribbed, RA
TPR024	Prophy Cup, Grey, 4 Web, Ribbed, RA
TPY000	Prophy Cup, Grey, Ribbed, RA
TPY002	Prophy Cup, Mint, 4 Web, Ribbed, RA
TPY002A	Prophy Cup, Grey, 4 Web, RA
TPY002E	Prophy Cup, Yellow, 4 Web, Ribbed,
TPY002G	Prophy Cup, Mint, 4 Web, RA
TPY003A	Prophy Cup, White, 6 Web, RA
TPY003C	Prophy Cup, Orange, 6 Web, RA
TPY809	Prophy Cup, White, Ribbed, SC
TPY805	Prophy Cup, Grey, Pointed, SC
PSO001	Prophy Cup, Pink, 6 Web, Snap-On
PSO002	Prophy Cup, Purple, 6 Web, Snap-On
PSO003	Prophy Cup, Blue, 6 Web, Snap-On
PSO004	Prophy Cup, White, 6 Web, Snap-On
PSO005	Prophy Cup, Yellow, 6 Web, Snap-On
PSO011	Prophy Cup, Pink, 4 Web, Ribbed, Snap-On
PSO012	Prophy Cup, Purple, 4 Web, Ribbed, Snap-On
PSO013	Prophy Cup, Blue, 4 Web, Ribbed, Snap-On
PSO015	Prophy Cup, Yellow, 4 Web, Ribbed, Snap-On
PSO019	Prophy Cup, Green, 4 Web, Ribbed, Snap-On
PSO020	Prophy Cup, Yellow, 4 Web, Ribbed, Snap-On
PSO021	Prophy Cup, Grey, 4 Web, Ribbed, Snap-On
PSO022	Prophy Cup, Green, 4 Web, Ribbed, Snap-On
PSO023	Prophy Cup, White, 4 Web, Ribbed, Snap-On
PSO024	Prophy Cup, Blue, 4 Web, Ribbed, Snap-On
PSO028	Prophy Cup, Blue, 4 Web, Ribbed, Snap-On
PSO029	Prophy Cup, Blue, 4 Web, Ribbed, Snap-On
EPR004A	Prophy Cup, White, 6 Web, Snap-On
EPR004B	Prophy Cup, White, 4 Web, Ribbed, Snap-On
TPR004	Prophy Cup, White, 4 Web, Snap-On
TPR004E	Prophy Cup, Blue, 4 Web, Ribbed, Snap-On
TPR004F	Prophy Cup, Yellow, 4 Web, Snap-On
TPR004K	Prophy Cup, Grey, 4 Web, Snap-On
TPR004L	Prophy Cup, Grey, 4 Web, Snap-On

### RECOMMENDED SPEEDS

- Maximum speed 10,000 rpm
- Keep within the speed range of 1,500 – 10,000 rpm for the best work results and to increase the service life of the instrument.

### INTENDED CLINICAL USE

The devices are indicated for use by qualified dental professionals for a wide range of restorative procedures, including the cleaning and polishing of natural teeth, gold and amalgam fillings, composites, compomers and glass ionomer cements for removing plaque and stains. The devices may be used in combination with prophylaxis polishing pastes which are themselves CE marked medical devices intended for use in the oral cavity. No other materials, substances or gases are used within the device.

Snap-On brushes and polishers are to be used with mandrel Ref: TPR010 or equivalent. Stoddard mandrels and brushes or polishers are proven to function together correctly by historical use with no adverse reports.

**Please refer to the Sterilisation and Instructions for Use for Single Use Devices for the correct procedures.**